

MiniVap™

Introducing the NEW evaporator for Porvair Sciences. The MiniVap™ replaces the Micro DS96 and joins the UltraVap in the state of the art 96 well sample evaporator units.



The MiniVap™ is designed for the research and development end of the market when individual plates need drying.

- Quicker dry down times than standard methods.
- Tests have shown that the MiniVap™ can evaporate 500µl of methanol in less than 10 minutes.
- Designed for any 96 well plate.
- Simple to install and operate.
- Easy adjustments of temperature, gas flow rates and needle depth into the wells.
- Compact foot print, will fit into all standard fume cupboards.

porvair
advanced materials
sciences

www.porvair-sciences.com

Unit 6
Shepperton Business park
Govett Avenue
Shepperton
Middlesex TW17 8BA
United Kingdom
Tel: Int +44 (0) 1932 240255
Fax: Int +44 (0) 1932 254393